



# Nanotechnology Enabled Image Guided Therapeutics in Lung Cancer

Kazuhiro Yasufuku MD, PhD

Director, Interventional Thoracic Surgery Program

Associate Professor of Surgery, University of Toronto

Division of Thoracic Surgery, Toronto General Hospital

*Personalizing Cancer Medicine in 2015*

*Feb 2<sup>nd</sup>, 2015*

---



# Disclosure

- Industry-sponsored grants
  - Educational and research grants from Olympus Medical Systems Corp.
- Consultant
  - Olympus America Inc.
  - Intuitive Surgical Inc.
  - Covidien
  - Johnson and Johnson
- Research Collaboration
  - Olympus Medical Systems Corp.
  - Novadaq Corp.
  - Veran Medical Technologies



## Lung Cancer

- Lung cancer remains the leading cause of cancer death in the Western world
- Early detection is key to improved survival
- The detection rate of early-stage lung cancer is anticipated to increase with the introduction of CT screening





## Lung Cancer – New Problems


- How can we efficiently localize small lung cancer during MIS?
- What are other minimally invasive therapeutic options for high risk pts with lung cancer?



2009



2011



2012






# Diagnostic Approach to Pulmonary Nodules


- Minimally Invasive Biopsy

- Bronchoscopic biopsy
- CT guided FNA



- Surgical biopsy

- VATS
- Thoracotomy





# Bronchoscopic Biopsy


- Transbronchial biopsy



- EBUS-GS



- Electromagnetic Guidance



- Virtual Navigation





## VATS (Video-assisted thoracoscopic surgery)

- Procedure of choice for surgical biopsy of peripheral pulmonary nodule
- Limitation
  - Identification of the nodule
  - Lack of digital palpation in small, non-solid deep nodules
  - May require conversion to thoracotomy



# Localizing Techniques - VATS





- Intraoperative imaging
  - CT
  - Thoracic ultrasound
- Preoperative CT guided marking
  - Liquid material (contrast media, colored adhesive agents, dyes)
  - Radionuclides
  - Wires (hookwires, microcoils)
- Preoperative bronchoscopic marking
  - Dye
  - Fiducials

*Radiology 2002; 225: 511-518*





# Wire Localization - Microcoil Localization




73/75 (97%) 4-24-mm nodules successfully removed  
low rate of intervention (3%) for procedural complications

*Radiology 2009; 250: 576-585*






# Microcoil Localization



1. CT guided micro coil placement



2. VATS microcoil detection



3. Fluoroscopy guided VATS wedge resection





4. Confirmation of microcoil



# Microcoil Localization – Toronto Experience

- First case October 2008
- 64 cases
- Complete resection with VATS in 62/64 cases (97%)
- 100% diagnostic yield






# GTx (Guided Therapeutics) Program





# GTx Surgery Overview



**TRIGOR = Translation Research Image Guided OR**



# TRIGOR A Capabilities



Cone-Beam CT



Dual Source- Dual Energy CT

MIS, Endoscopic Technology






# TRIGOR A - GTx OR



# Multi-Modality Surgical Guidance

## Pre-Operative Imaging



## Intraoperative CBCT




## Surgical Planning

## Optical Imaging


## Surgical Tool Tracking




## Intraoperative localization



## Image assistance during MIS/Robotic Surgery



## Real time monitoring of minimally invasive thoracic intervention





# Optical Surgical Navigation

1. Preoperative CT imaging



2. Cone-beam CT nodule localization




3. Image registration with surgical navigation


Pre-Operative Imaging



Intraoperative CBCT



Surgical Tool Tracking





## GTx OR – Image guided Transbronchial Interventions






## GTx OR – VATS localization







# ICG - SPY Localization








# Image-guided Localization Platform for Minimally Invasive Lung Cancer Surgery





# Multi-modal liposomal nanoparticle(C800) - ICG Liposome



Size (nm)	Iodine concentration (mg/mL)	ICG concentration (µg/mL)
99.2 ± 1.8	53.2 ± 1.4	94.6 ± 7.8

The phospholipid nanoparticle, coated with polyethylene glycol, encapsulates ICG and CT contrast (iohexol, labelled IOX). The prolonged intravascular half-life allows for longitudinal CT and NIR imaging.






## Image Assistance – 3D Reconstruction

Pre-injection

5 min post-injection

3 days post-injection




Administered CF800 continuously remains in both vasculature and the tumor in a rabbit, allowing for successful 3D reconstruction even at 3 days post injection







## NIR imaging of Lung Cancer



4 days post-injection of CF800





# Ultra-minimally Invasive multi-modal image guided photothermal ablation of lung cancer





# Porphysomes : Liposomes like bilayer porphyrin-phospholipid



ARTICLES

PUBLISHED ONLINE: 20 MARCH 2011 | DOI:10.1038/NMAT2986

nature  
materials

## Porphysome nanovesicles generated by porphyrin bilayers for use as multimodal biophotonic contrast agents



Jonathan F. Lovell<sup>1,2</sup>, Cheng S. Jin<sup>1,2</sup>, Elizabeth Huynh<sup>2,3</sup>, Honglin Jin<sup>2,3</sup>, Chulhong Kim<sup>4</sup>, John L. Rubinstein<sup>3,5</sup>, Warren C. W. Chan<sup>1</sup>, Weiguo Cao<sup>6</sup>, Lihong V. Wang<sup>4</sup> and Gang Zheng<sup>1,2,3\*</sup>




Works as both fluorophore and photo-enhancer




# EPR effect; Enhanced Permeability and Retention Effect



Cancer vasculature




Normal vasculature



- 1, Adv Drug Deliv Rev, 2011
2. Drug Dis Today, 2006





## Porphysome distribution in Orthotopic lung cancer Xenograft models






## Porphysome thermal effect on resected VX2 tumor tissue ( ex-vivo)





# Development of prototype bronchoscope for visualization of porphyrin fluorescence




## Emission (bronchoscope side) :

A long pass filter : 678-1000nm

## Excitation (light source side):

Red diode of 10 mW output  
A band-pass filter : 650-670nm






# In-vivo model

## Prototype fluorescence bronchoscope


WL




0 time



24hr




48hr








# Experimental Design




## Tumor-growth






# Porphysome enhanced transbronchial PTT

Laser Control




H&E staining




NADH staining


•Ablated area: < 2 mm diameter





Lung Slide



H&E Staining



NADH Staining




Porphysome-PTT





## Ultra-minimally invasive multi-modal image guided photothermal ablation of lung cancer





# Summary

- Advances in nanotechnology and image guidance will enable intraoperative localization of small peripheral nodule and also assist surgeons during MIS
- New transbronchial ablation technologies are in development and can potentially be used for minimally invasive treatment of early stage lung cancer





# Acknowledgement

## Toronto General Hospital

Doctors of Thoracic Surgery (DOTS)

Shaf Keshavjee

Tom Waddell

Gail Darling

Andrew Pierre

Marc de Perrot

Marcelo Cypel



## Latner Thoracic Surgery Laboratory

Yasufuku Laboratory

Hironobu Wada

Tatsuya Kato

Priya Patel

Spencer Hu

Daiyoon Lee

Judy McConnell

Alexandria Grindlay

Muhammed Hassan





---

# Acknowledgement

## **Department of Medical Biophysics**

Gang Zheng  
Brian Wilson  
Warren Chan  
Jonathan Lovell

## **Faculty of Medicine**

David Hwang  
Ren-Ke Li

## **TECHNA Institute**

David Jaffray  
Robert Weersink  
Jonathan Irish  
Harley Chan  
Michael Dunne  
Jimmy Qiu

## **Chiba University**

Department of Thoracic Surgery  
Prof Yoshino

## **Kochi University**

Department of Surgery  
Prof Orihashi

## **Hokkaido University**

Department of Cardiovascular and  
Thoracic Surgery  
Prof Matsui  
Prof Kaga





Division of Thoracic Surgery  
Toronto General Hospital  
University Health Network

Kazuhiro Yasufuku, MD, PhD, FCCP  
[kazuhiro.yasufuku@uhn.ca](mailto:kazuhiro.yasufuku@uhn.ca)

Thank you